

GEROtherm® UNI-FIX

As a weight connector for universal weight attachment and the fastening of double U GEROtherm® probe feet 32–43 mm.

Max. load: 100 kg.

No mechanical safeguard

Length approx. 230 mm

Art. No.: 104606

Quantity: Units

UNI-FIX 2020

Length approx. 245 mm

Art. No.: 113284

Quantity: Units

GEROtherm® weight connector de50mm

As a weight connector for universal weight attachment and the fastening of double U GEROtherm® probe feet 50 mm. Matching with additional weight 20 kg (104644)

Max. load: 300 kg.

Length approx. 340 mm

Art. No.: 116841

Quantity: Units

GEROtherm® PUSH-FIX

Impact-resistant sleeve and weight connector for the installation of GEROtherm® probe feet (32–43 mm) with linkage support.

Max. impact: 2,000 kg.

Single U 32 mm

Length approx. 370 mm

Outside diameter with deflection: approx. 94 mm

Art. No.: 108139

Quantity: Units

Single U 35 mm REX

Length approx. 370 mm

Outside diameter with deflection: approx. 96 mm

Art. No.: 108140

Quantity: Units

Single U 40 mm

Length approx. 385 mm

Outside diameter with deflection: approx. 112 mm

Art. No.: 108141

Quantity: Units

Single U 43 mm REX

Length approx. 385 mm

Outside diameter with deflection: approx. 112 mm

Art. No.: 108143

Quantity: Units

Double U 32 mm
Length approx. 370 mm
Outside diameter with deflection: approx. 100 mm
Art. No.: 108144
Quantity: Units

Double U 35 mm REX
Length approx. 370 mm
Outside diameter with deflection: approx. 105 mm
Art. No.: 108146
Quantity: Units

Double U 40 mm
Length approx. 385 mm
Outside diameter with deflection: approx. 118 mm
Art. No.: 108147
Quantity: Units

Double U 43 mm REX
Length approx. 385 mm
Outside diameter with deflection: approx. 120 mm
Art. No.: 108149
Quantity: Units

GERO[®]therm[®] INTERFACE

for simple incorporation of the impact-resistant sleeve in combination with the GERO[®]therm[®] PUSH-FIX
Length approx. 635 mm
Art. No.: 108157
Quantity: Units